

NeuroTherm®

Curved Disposable RF Cannula

Always
a step
ahead...



Curved Sharp or Blunt

Our curved cannula are available in sharp and blunt tips, in a variety of lengths, active tip sizes and gauges ensuring we have the right one for you

Compatibility

Our straight and curved cannula are safely compatible with all major brands of RF generators

Customization


Still can't find what you are looking for? We can manufacture custom lengths, active tip sizes and gauges to meet your specific needs

NeuroTherm®
30 Upton Drive, Suite 2
Wilmington, MA 01887
USA
Tel: +1 888 655 3500
Fax: +1 978 658 2378



NeuroTherm®
429 Brighton Road, Croydon
Surrey, CR2 6EU
United Kingdom
Tel: +44 20 8660 4374
Fax: +44 20 8660 9417

NeuroTherm®
Modemstraat 20
1033 RW Amsterdam
The Netherlands
Tel: +31 204 702 702
Fax: +31 204 702 701



www.neurotherm.com

NeuroTherm®

The Disposable Curved RF Cannula

Our disposable curved sharp and blunt RF cannula are designed for use with radiofrequency electrodes and generators in (P)RF procedures. Our cannula have machine electro-polished tips for easy introduction. The curve allows for better cannula steerability with the red dot on the ergonomic hub indicating curve direction at the tip.



Length	Tip	Gauge	Part #
5 cm Curved Sharp			
5 cm	2 mm	22	C-502
5 cm	5 mm	22	C-505
5 cm	5 mm	20	C-505-20
5 cm	5 mm	18	C-505-18
5 cm	10 mm	22	C-510
5 cm	10 mm	20	C-510-20
5 cm	10 mm	18	C-510-18
10 cm Curved Sharp - SMK			
10 cm	5 mm	22	C-1005-S
10 cm	5 mm	20	C-1005-S-20
10 cm	5 mm	18	C-1005-S-18
10 cm	10 mm	22	C-1010-S
10 cm	10 mm	20	C-1010-S-20
10 cm	10 mm	18	C-1010-S-18
15 cm Curved Sharp - SMK			
15 cm	5 mm	20	C-1505-S
15 cm	5 mm	18	C-1505-S-18
15 cm	10 mm	20	C-1510-S
15 cm	10 mm	18	C-1510-S-18
10 cm Curved Sharp - RFK			
10 cm	5 mm	22	C-1005-R
10 cm	5 mm	20	C-1005-R-20
10 cm	5 mm	18	C-1005-R-18
10 cm	10 mm	22	C-1010-R
10 cm	10 mm	20	C-1010-R-20
10 cm	10 mm	18	C-1010-R-18
15 cm Curved Sharp - RFK			
15 cm	5 mm	22	C-1505-R
15 cm	5 mm	18	C-1505-R-18
15 cm	10 mm	20	C-1510-R
15 cm	10 mm	18	C-1510-R-18
5 cm / 10 cm Curved Blunt			
5 cm	5 mm	22	C-505-B-22
10 cm	10 mm	20	C-1010-B
10 cm Curved Blunt - SMK			
10 cm	5 mm	22	C-1005-S-B-22
10 cm	5 mm	20	C-1005-B-S
10 cm Curved Blunt - RFK			
10 cm	10 mm	22	C-1005-R-B-22
10 cm	10 mm	20	C-1005-R-B
15 cm Curved Blunt - SMK			
15 cm	5 mm	20	C-1505-S-B
15 cm	10 mm	20	C-1510-S-B
15 cm Curved Blunt - RFK			
15 cm	10 mm	20	C-1505-R-B
15 cm	10 mm	20	C-1510-R-B

Distributed by:

www.neurotherm.com



ISO 13485 Certified

NeuroTherm is a registered trademark of NeuroTherm Inc.

SS-108 Rev. 1 06/10