

Always
a step
ahead...

Spinal Needles

Expanding your choices for pain management

Curved Tip Needles

Blunt (For Passing through Sensitive Tissue)

Part Number	Length	Gauge
B-420	4	20
B-422	4	22
B-620	6	20
B-622	6	22

Curved Quincke Spinal Needles

Sharp (Curved for Steering)

Part Number	Length	Gauge
C-222	2	22
C-420	4	20
C-422	4	22
C-620	6	20
C-622	6	22

Straight Spinal Needles

Part Number	Length	Gauge
S-422	4	22
S-522	5	22
S-622	6	22
S-722	7	22
S-822	8	22



NeuroTherm®
30 Upton Drive, Suite 2
Wilmington, MA 01887
USA
Tel: +1 888 655 3500
Fax: +1 978 658 2378



NeuroTherm®
429 Brighton Road, Croydon
Surrey, CR2 6EU
United Kingdom
Tel: +44 20 8660 4374
Fax: +44 20 8660 9417

NeuroTherm®
Modemstraat 20
1033 RW Amsterdam
The Netherlands
Tel: +31 204 702 702
Fax: +31 204 702 701

Disposable Discography Kits

The NeuroTherm® Discography Kit,

Consists of two needles. The needle hub design is based on the SMK RF Injection needle, and has a luer lock to facilitate the fitting of a syringe

The outer needle (20g/18g) with stylet is inserted under fluoroscopy until it is in the annulus fibrosus.

The stylet is removed and the second needle (25g/22g) is inserted into the outer needle and advanced into the nucleus pulposus.

A contrast medium is injected into the disc, which will demonstrate if tears, scarring and/or changes in the lining of the disc are present. Additionally, disc pain may be provoked by the injection of the contrast medium (Provocative Discography), which is a second indicator that a disc is the source of the pain.

Part Numbers:

Part Number	Description
DK-1822	Discography Kit - 18g x 3.5" introducer x 22g x 6.0" chiba needles
DK-1822-80	Discography Kit - 18g x 3.5" introducer x 22g x 8.0" chiba needles
DK-1822-68	Discography Kit - 18g x 6.0" introducer x 22g x 8.0" chiba needles
DK-2025	Discography Kit - 20g x 3.5" introducer x 25g x 6.0" chiba needles
DK-2025-8	Discography Kit - 20g x 3.5" introducer x 25g x 8.0" chiba needles

NeuroTherm® combines a strong reputation for innovative designs that are safe, reliable and easy to use, with a remarkable level of service, so physicians can treat their patients more safely and effectively.

Contact us today, to learn more about our products and services.

Distributed by:

www.neurotherm.com



ISO 13485 Certified

0086
NeuroTherm is a registered trademark of NeuroTherm Inc.