

Straight Disposable RF Cannula

Always
a step
ahead...



Best Selection

The widest variety of cannula lengths, active tip sizes and gauges ensures we'll always have the right cannula available for you

Compatibility

Our straight and curved cannulae are safely compatible with all major brands of RF generators

Customization


Still can't find what you are looking for? We can manufacture custom lengths, active tip sizes and gauges to meet your specific needs

NeuroTherm[®]
30 Upton Drive, Suite 2
Wilmington, MA 01887
USA
Tel: +1 888 655 3500
Fax: +1 978 658 2378



NeuroTherm[®]
429 Brighton Road, Croydon
Surrey, CR2 6EU
United Kingdom
Tel: +44 20 8660 4374
Fax: +44 20 8660 9417

NeuroTherm[®]
Modemstraat 20
1033 RW Amsterdam
The Netherlands
Tel: +31 204 702 702
Fax: +31 204 702 701



The Disposable Straight RF Cannula

Our disposable straight RF cannula is designed for use with radiofrequency electrodes and generators in (P)RF procedures. Our cannulae have machine sharpened, electro-polished tips for easy introduction. The ergonomic hub allows for easy advancement and the quality insulation prevents RF leakage.



5 cm Straight RF Cannulae			
Length	Tip	Gauge	Part #
5 cm	2 mm	22	S-502
5 cm	2 mm	21	S-502-21
5 cm	5 mm	22	S-505
5 cm	10 mm	22	S-510
10 cm Straight RF Cannulae			
10 cm	2 mm	22	S-1002
10 cm	5 mm	22	S-1005
10 cm	5 mm	21	S-1005-21
10 cm	5 mm	20	S-1005-20
10 cm	5 mm	18	S-1005-18
10 cm	10 mm	22	S-1010
10 cm	10 mm	21	S-1010-21
10 cm	10 mm	20	S-1010-20
10 cm	10 mm	18	S-1010-18
15 cm Straight RF Cannulae			
15 cm	2 mm	20	S-1502
15 cm	5 mm	20	S-1505
15 cm	5 mm	18	S-1505-18
15 cm	10 mm	20	S-1510
15 cm	10 mm	18	S-1510-18
20 cm Straight RF Cannulae			
20 cm	5 mm	20	S-2005
20 cm	10 mm	20	S-2010
* For use with Smith & Nephew™			
10 cm	5 mm	22	S-1005-SN
10 cm	5 mm	20	S-1005-SN-20
10 cm	10 mm	22	S-1010-SN
10 cm	10 mm	20	S-1010-SN-20
15 cm	5 mm	20	S-1505-SN
15 cm	10 mm	20	S-1510-SN

Distributed by:

www.neurotherm.com



ISO 13485 Certified